



Integration of continuous upstream and downstream operations in mAb production

Intellectual Property Notice: The Biopharma business of GE Healthcare was acquired by Danaher on 31 March 2020 and now operates under the Cytiva™ brand. Certain collateral materials (such as application notes, scientific posters, and white papers) were created prior to the Danaher acquisition and contain various GE owned trademarks and font designs. In order to maintain the familiarity of those materials for long-serving customers and to preserve the integrity of those scientific documents, those GE owned trademarks and font designs remain in place, it being specifically acknowledged by Danaher and the Cytiva business that GE owns such GE trademarks and font designs.

cytiva.com

GE and the GE Monogram are trademarks of General Electric Company. Other trademarks listed as being owned by General Electric Company contained in materials that pre-date the Danaher acquisition and relate to products within Cytiva's portfolio are now trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate. All other third-party trademarks are the property of their respective owners.

© 2020 Cytiva

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit [cytiva.com/contact](https://www.cytiva.com/contact)

CY14323-10Jun20-PT

Integration of continuous upstream and downstream operations in mAb production

Andreas Castan, Thomas Falkman, Eric Fäldt, and Annika Forss

GE Healthcare Bio-Sciences AB, Björkgatan 30, 751 84 Uppsala, Sweden

Introduction

Process intensification is gaining interest as a strategy to reduce production costs, while improving product quality and throughput in the manufacturing of biopharmaceuticals. For a competitive production process, continuous or semi-continuous upstream and downstream processing can be employed. Compared with a process performed in batch runs, continuous processing allows for increased capacity utilization and eliminates or minimizes the need for intermediate hold-up steps. We describe the integration of a high-performing upstream cell culture process with downstream

purification utilizing emerging technologies such as periodic counter-current (PCC) chromatography and straight-through processing (STP).

In this case study, we demonstrate the feasibility of integrated upstream and downstream mAb processing performed in a continuous manner. The developed process shows a performance equivalent to traditional processing performed in batch runs.

Results

Upstream: perfusion medium development

A high-performing cell culture medium supporting a low cell-specific perfusion rate (CSPR) is the key to developing a sustainable perfusion process. To develop a high-performing perfusion medium, we studied combinations in an existing medium platform, consisting of medium and feeds (Fig 1).




Fig 1. Study design for developing a high-performing perfusion medium from an existing fed-batch platform.

High-performing perfusion media can be developed from an existing medium platform, the perfusion rate can be used as a tuning fork for product quality (Fig 2 and Table 1).




Fig 2. Confirmation of perfusion media formulation from the steady-state approach at 50 MVC/mL and 1 RV/d (= CSPR 20).

Table 1. Quality attributes impacted by the perfusion rate (quality attributes not impacted were: glycosylation profile and TNF- α binding)

| Analytical technology | Analyte | Fed-batch | 20 pL/c/d | 43 pL/c/d | 77 pL/c/d | 90 pL/c/d |
|-----------------------|-------------------|-----------|-----------|-----------|-----------|-----------|
| CIEX | Acidic variants | > 60% | 25% | 23% | 19% | 12% |
| CIEX | Alkaline variants | 3% | 2% | 4% | 4% | 4% |
| SEC | Aggregate | 1% | 0.4% | 0.3% | 0.2% | 0.5% |

Downstream: capture by 3C PCC

The capture step was performed running on a 3C PCC setup using HiScreen™ columns prepacked with MabSelect SuRe™ chromatography resin. The principle is shown in Figure 3. The performance and robustness was evaluated, running 10 cycles, with respect to purity and yield. Loading was controlled by dynamic control, which enables system operation under process conditions where either feed concentration or chromatography medium capacity varies (Fig 4).




Fig 3. (A) Working principle of 3C-PCC; (B) schematic view of UV curves during loading.




Fig 4. Demonstration of the dynamic control. The dynamic control was challenged by performing a feed gradient run from 1 to 2 g/L IgG. The breakthrough curves follow the gradient change closely, the amount of product eluted from each column remained constant (data not shown).

Robust performance and minimal variation in amount mAb, aggregate, and HCP between columns was determined (Fig 5). Results were comparable to results from batch protein A capture. A 56% increase in capacity compared with batch mode was observed (67 g/L, compared with 43 g/L, loaded onto each column).




Fig 5. The elution pool from each individual column and cycle was analyzed with respect to the amount mAb, aggregate content, and host cell protein. Results from each column in the graphs are presented according to column number.

Downstream: polishing by STP

A polishing STP step was developed including the following steps:

- Individual studies on the different chromatography resins
- Conditioning study based on results from studies on Capto™ S ImpAct and Capto adhere resins
- Verification runs with Capto S ImpAct and Capto adhere connected in series with in-line conditioning between the steps

Figure 6 illustrates the set-up used. Table 2 shows the quality indicating analytical data.




Fig 6. Process options investigated during verification runs.

Table 2. Results from the STP step: B/E-FT and FT-FT process alternatives. The starting material contained 2% aggregates, 1277 ppm HCP, and 3 ppm leached ligand.

| Process step | Yield over two steps (%) | Aggregates (%) | HCP (ppm) | Leached ligand (ppm) |
|--|--------------------------|----------------|-----------|----------------------|
| Capto S ImpAct (B/E)- Capto adhere (FT) | 89 | 0.8 | 16 | < 1 |
| Capto S ImpAct (B/E)- Capto adhere (FT) | 90 | 0.7 | 18 | < 1 |
| Capto S ImpAct (FT)- Capto adhere (FT) | 88 | 0.7 | 40 | < 1 |

Yield as well as amount aggregate, HCP, and leached ligand were comparable with results from batch runs.

Abbreviations

| | |
|--------|--|
| 3C PCC | three-column periodic counter-current chromatography |
| B/E | bind-and-elute mode |
| CIEX | cation exchange chromatography |
| CSPR | cell-specific perfusion rate |
| DoE | design of experiments |
| FT | flow-through mode |
| HCP | host cell protein |
| PAT | process analytical technology |
| PCC | periodic counter-current chromatography |
| RV | reactor volume |
| SEC | size exclusion chromatography (gel filtration) |
| STP | straight-through processing |

Conclusions

- The cell culture medium is critical for the development of high-performing perfusion processes.
- High-performing perfusion media can be developed from existing fed-batch platforms, allowing CSPR's of 20 pL/c/d.
- The perfusion rate can be used as a tuning fork for product quality.
- Continuous protein A capture increased productivity compared with a batch process and delivered robust performance.
- The dynamic control is an example of feed-back control and integrated PAT.
- An STP step with two columns connected with in-line adjustment in-between showed similar yield and purity as batch runs.